



Patient Management



Import



Export



advanced G.I. diagnostic software

High Resolution Impedance
Manometry (HRiM[®])

Swallow Atlas



800-558-6408
303-470-7020




www.sandhillsci.com




clinicalsupport@sandhillsci.com
technicalsupport@sandhillsci.com

SANDHILL
SCIENTIFIC
a medtronic company

Navigation icons: Home, Refresh, Save, Print, Zoom (00:30), Full Screen, Split Screen, Window Management, Help.



01:16 01:18 01:20 01:22 01:24 01:26 01:28 01:30 01:32 01:34 01:36 01:38 01:40 01:42 01:44







Normal Swallow





Navigation toolbar with icons for undo, redo, save, print, zoom (00:30), and various view controls.



Ready







Pressure

PHARYNX 150

UES 9.9 cm

BODY 21.3 cm

LES 3.9 cm CD 5.5 cm


mmHg


Impedance


5000


ohms

05:34 05:36 05:38 05:40 05:42 05:44 05:46 05:48 05:50 05:52 05:54 05:56 05:58 06:00 06:02











Navigation icons: Home, Refresh, Save, Print, Zoom (00:30), and various view controls.





PHARYNX
UES
BODY 24.8 cm
LES 3.8 cm
CD

Pressure 150 mmHg
Impedance 5000 ohms

12:38 12:40 12:42 12:44 12:46 12:48 12:50 12:52 12:54 12:56 12:58 13:00 13:02 13:04 13:06









Pressure

150
30
20
-10
mmHg

Impedance


5000
650
0
ohms

PHARYNX
UES
5.7 cm
BODY
22.3 cm
LES
CD



Sandhill University—training and education at a higher level.

Providing great products is just the first step in a long term relationship with our customers...we believe that training and education are the key components to even the best equipment, enabling you to provide outstanding patient care. That's why at Sandhill University, we have developed a comprehensive set of options to meet your needs, including: Sandhill University/Denver—Sandhill University/Online—Webinars—Cybercoaching.



Manometry Testing



Reflux Monitoring



Esophageal Dilators



CleanFreak® Cleaning Products



Endoscope Care & Accessories

Sandhill Scientific—a global leader in G.I. diagnostics.

Since 1981, Sandhill Scientific has been recognized as a global leader in the field of G.I. Diagnostics. With a long and rich history in the G.I. space, we have produced many of the innovations used today to help enhance the diagnostic yield of reflux and manometry.

With our Technical Research & Training Center in Denver, Sandhill will continue to innovate, blending science and user feedback to develop products and technologies for you, our customers. See the complete line of GI solutions at **SandhillSci.com**.



Corporate Headquarters

102 East Keefe Avenue
Milwaukee, WI 53212 USA
800.558.6408 or 414.265.7635
Fax 414.265.7628

sales@sandhillsci.com
orders@sandhillsci.com

Technical Research & Training Center

9150 Commerce Center Circle #500
Highlands Ranch, CO 80129 USA
800.558.6408 or 303.470.7020
Fax 303.470.2975

technicalsupport@sandhillsci.com
clinicalsupport@sandhillsci.com
clinicaleducation@sandhillsci.com